



TELFA[™] AMD[™]

ANTIMICROBIAL DRESSINGS

- **Low Cost** — Pennies per day
- **Convenient** — Works in moist or dry conditions
- **Simple** — Compliant with existing protocol



ZONE OF INHIBITION



- **TELFA AMD Dressing vs. untreated dressing against Staphylococcus aureus.***
- **TELFA AMD Dressing vs. untreated dressing against Escherichia coli.***
- **TELFA AMD Dressing vs. untreated dressing against Candida albicans.***

* Dressing in contact with seeded agar plate modulated at 10^6 CFU/ml.

Kendall AMD Dressings



Did you know?

500,000 SSI's occur annually in the USA*?

Estimated average extra cost due to SSI is \$3,500 per patient*?

Nosocomial Infections are the fourth leading cause of Hospital deaths*?


- **TELFA is the preferred non-adherent dressing for surgical sites.**
- **Protect your surgical sites with TELFA AMD.**

* According to CDC

AMD is effective against:

Staphylococcus aureus
Staphylococcus epidermidis
Pseudomonas aeruginosa
Escherichia coli
Candida albicans
Staphylococcus coagulase
Proteus mirabilis
Serratia marcescens
Enterobacter cloacae
Klebsiella pneumoniae
Enterococcus faecalis

Graph: AMD Efficacy



*0/24/48 at 10⁶ inoculation levels for challenge organisms

Ordering Information

TELFA AMD Pads

Code	Description	Per Tray	Ship Case
7662	3" x 4" Pad, Sterile 1's	50	900
7663	3" x 8" Pad, Sterile 1's	50	600

TELFA AMD Island Dressings

Code	Description	Per Tray	Ship Case
7665	4" x 5" Island, Sterile 1's	25	200
7666	4" x 8" Island, Sterile 1's	25	100
7667	4" x 10" Island, Sterile 1's	25	100
7668	4" x 14" Island, Sterile 1's	25	50

™ Trademarks of Tyco Healthcare Group LP

tyco
Healthcare

Kendall

15 Hampshire Street
Mansfield, MA 02048
1-800-962-9888
508-261-8000

www.kendallhq.com

For International locations please see
www.tycohealthcare.com

H-5656 10M ©Tyco 9/2004